

SPORT INJURIES

BIOPTRON LIGHT THERAPY

Duration: 4-6 minutes

Tennis Elbow, golfers elbow Tendonitis	R300 for 3 sessions
Calf and hamstring injuries Back pain	or
Swelling Spasms and knots Neck pain	R500 for 6 sessions
Pre and post training stiffness Sprains, strains Soft tissue injuries of muscles	

PAIN TREATMENT

BIOPTRON LIGHT THERAPY

Duration: 6-10 minutes

Rheumatology: Osteoarthritis, rheumatoid arthritis	R300 for 3 sessions
Physiotherapy: Low back pain, shoulder and neck pain	or
	R500 for 6 sessions

WOUND HEALING

BIOPTRON LIGHT THERAPY

Burns	R300 for 3 sessions
Wound healing after trauma (injuries)	or
Grafting	R500 for 6 sessions

ANTI AGING

BIOPTRON LIGHT THERAPY

Duration: 10 minutes

Stimulate Skin Renewal Collagen Production	R300 for 3 sessions
Reduce Appearance of Pores Skin Smoothness	or
Enhance Radiance Diminish Fine Lines	R500 for 6 sessions
Improves Firmness and Elasticity Microcirculation Enhancement Cellular Energy Support	

DERMATOLOGICAL DISORDERS AND SKIN PROBLEMS

BIOPTRON LIGHT THERAPY

Duration: 4-6 minutes

Psoriasis Superficial bacterial infections	R300 for 3 sessions
Acne (rosacea) Atopic dermatitis	or
	R500 for 6 sessions



SLIMMING

Pressotherapy

Stimulate lymphatic system	R200
Improves circulation	for 25
Promotes detoxification	minutes
Reduces Swelling Alleviates water retention Enhances body contouring	

Far Infrared Therapy

Stimulates metabolism	R200
Improves blood circulation	for 25
Burns fat and boosts metabolism	minutes
Enhances nutrient absorption Increases collagen production	

EMS (Electrical Muscle Stimulation)

Stimulate muscle contraction .Tones and strengthens muscle	R200
Body shaping	for 25
Relieves muscle tension	minutes
Aids post-workout recovery	